



Class 3 Invalid Carriage

- Is an electric powerchair or scooter.
- Can be driven along roads and can be used on dual carriage ways with speeds of under 50mph.
- Usually has a top speed of 8pmh.
- Can be used on pavements but must be restricted to 4mph when they do so.
- Has a minimum rider age of 14 years old.
- A class 3 Road scooter MUST be registered with the DVLA.

Class 3 mobility scooters are subject to more stringent rules than a pavement scooter. To be legally be driven on the road a class 3 scooter must also have the following:

- A maximum unladen weight of 150kg (or a maximum unladen weight of 200kg when carrying necessary user equipment, for example, medical supplies).
- A maximum width of 0.85 metres.
- A device to limit its speed to 4mph.
- A maximum speed of 8mph.
- An efficient braking system.
- Front and rear lights and reflectors.
- Direction indicators able to operate as a hazard warning signal.
- An audible horn.
- A rear view mirror.
- An amber flashing light if it's used on a dual carriageway.